




LIBERO PDF Logger for Cryogenic Container Monitoring

we prove it.

LIBERO Cryogenic Monitoring Solution



LIBERO Cryogenic Monitoring Kit




Suitable for most
cryogenic
container types




LIBERO Te1-PY
well-protected in
stainless steel bracket

- Data logger with stainless steel bracket is mounted on top of the lid; temperature probe is mounted through the cork into the upper part of the vapor phase area.
- Send us the lid (cork) or the complete container in order to mount the bracket and the probe for you.

Principle of Mounting on a Cryo Shipper



Data logger and Bracket



LIBERO PDF Logger Te1-P for external Pt100 probes

Data logger can be plugged into any USB port of a PC, MAC or Linux to generate a evaluation report, no software required

Part No	Description	Log interval min.	Battery
800021	LIBERO Te1-PS, single use	1 minute	100 days
800020	LIBERO Te1-PY, multiple use	1 minute	400 days
800019	LIBERO Te1-P, multiple use	5 minutes	2 years

Alarm Bracket for Pt100 Probe

This is mainly used for fixed installations. The alarm bracket can generate visual and audible alarm. It can trigger external alarms.

Part No	Description
800915	Alarm Bracket

Stainless steel bracket for Libero and Pt100 Probe

Mostly used for mobile cryo shippers, lockable by padlock
Data logger is well protected. Sensor is connected in a terminal.

Part No	Description
800911	Stainless Steel Bracket with terminals
800914	Stainless Steel Bracket to protect the Alarm Bracket


Z- Mounting Plate

To fix the bracket 4402-A to the side of a cork

Part No	Description
800912	Z-Mounting Plate

Overview CRYO Temperature Probes

-200°C .. +200°C



ELPRO CRYO Temperature Probes


All types are available in different lengths and with calibration certificate at -196°C

Type	Diameter
PTFE / Stainless Steel	5.0mm
PTFE / Stainless Steel	4.4mm
Stainless Steel, bendable	3.0mm
Stainless Steel, bendable	2.0mm
PTFE	2.2mm

Cryo Temperature Probes -200°C .. +200°C

Cryo Probe ML Type, Stainless Steel -200°C .. +200°C, bendable

(Example shown with 2 bends, delivered straight)



ELPRO CRYO Probe ML (Stainless Steel) D=3mm
bendable for easy adjustment to application
Pt100, Class A, temperature range: -200°C .. +200°C

Part No.	Text	Length L	Cable C
800685	Cryo Probe ML D=3mm	300mm	50mm
800686	Cryo Probe ML D=3mm	370mm	100mm
800687	Cryo Probe ML D=3mm	200mm	50mm
800688	Cryo Probe ML D=3mm	450mm	150mm

Additional Types of Cryo Probes are available on request

- With M12 connector and extension cable
- Made with PTFE cable in different diameters



LIBERO PDF Logger for external Probe and Alarm

- Automatic built-in PDF Generator
- Single & Multiple use versions (400 days/2 years)
- Memory: 16'000 temperature data points
- Measuring range: -200°C ..+200°C
- High accuracy: +/- 0.2° (depending on range)
- Large, multifunctional LCD display
- Alarm: 2 limit values or 5 alarm/time ranges
- Alarm indicator and external Alarm Output
- Function Keys: TRANSIT and ARRIVED
- Detailed evaluation of embedded data possible with free elproVIEWER software



Libero PDF Report No 89568811 (Asia_AB02 20070803194130 10000624.pdf)

Test Switzerland_CCM Report



Additional Information

Shipping No:

IMPORTANT!

A) EMAIL this LIBERO PDF-File immediately to: info@elpro.com
B) Send Logger back to: Elpro-Buchs AG, CH-8470 Buchs, Switzerland


Device Configuration

Type:	Libero Ti1 V1.15	Inspection Range:	Last Transit / Arrived
Logger ID:	10000624	Current State:	Logging Transit
Log Interval / Duration:	3 m / 33.3 d	Remaining Battery:	372 d
Log Mode:	Loop	Logger Start:	03.07.2007 16:42:18
Report Time Base:	GMT +02:00		
Configured by:	ELPRO067/ablschhof, 03.07.2007 16:41:03		

Alarm Zones	Allowed Time	Total Time	No. Violations	Status
Z1: over 35.0 °C	15 m (min)	0 s	0 / unlim.	OK
Z2: over 25.0 °C	3.0 d (cum)	6.8 h	3 / unlim.	OK
Z3: over 8.0 °C	5.0 d (cum)	34.8 h	2 / unlim.	OK
Z4: 2.0 to 8.0 °C	unlimited	36.1 h		
Z5: below 2.0 °C	24.0 h (cum)	49.3 h	1 / unlim.	ALARM
Z6: below 0.0 °C	5 m (min)	49.3 h	1 / unlim.	ALARM

Logging Results


Highest Temperature:	31.7 °C	Transit Start at:	31.07.2007 07:30:41
Lowest Temperature:	-26.7 °C	Arrived at:	not available
Average Temperature:	-3.5 °C	File created:	03.08.2007 19:41:30




Report created by Libero, the PDF Data Logger

Just plug LIBERO into any USB port:
PDF/A report, automatically created

Mounted on lid of MVE CryoShipper, -XC and -Mini




ww.elpro.com
AB June 2010



Mounting for smaller Lid Diameter, Bracket fixed to Tank, Vapour shipper: **MVE Doble 11**

Principle of mounting
Libero Cryo Monitor on MVE Doble 11
Cryoshipper



Elpro-Buchs AG, Switzerland
June 2010, AS
www.elpro.com





Mounting for Taylor Wharton Cryo Express CX 100

Libero Cryo Monitoring TW CX100



Mounting with fixation plate on side

Measurements
all in [mm];

www.elpro.com
AB June 2010




Taylor-Wharton / Pangas CX100-B-11M
Length of probe: $L = 38 + 136 + 80 + 27 + 39$
= 320 (incl. Endpart) Selection of Probe: 3148-B



The complete LIBERO PDF Logger is fixed to the cork and still fits into the standard transportation case .

Mounting on AIR LIQUIDE Voyageur 5





LIBERO Alarming Bracket for fixed Cryo Containers



Cryo Container is equipped with a temperature probe going inside the LN2 tank. When the container is back at its place in the laboratory, it can be connected to the LIBERO Alarm Bracket. This device detects the Alarms generated by the LIBERO PDF Logger and starts beeping. With the internal dry contact, a telephone auto dialler may be connected or any connection to the building alarming system.


Temperature Sensor goes into Cryo Vessel



Handling and System Components

At departure:


1. Check LIBERO display:
If it is already in «TRANSIT», press the «STOP/Arrived» button for 3 seconds.
2. **To start a new transport section,**
press «START/Transit» button for 3 seconds:



Logging status changes to «TRANSIT» and temperature monitoring becomes active.


At destination:

1. **To mark the end of a transport section,** press the black «STOP/Arrived» button for about 3 seconds:



Logging status changes to «ARRIVED».

2. Disconnect LIBERO from Terminal and plug LIBERO into the USB port of your PC.
3. Within few seconds, the evaluation report (PDF/A file) will automatically be generated and presented in the explorer window.



Costs and required Material

is depending on the type of cryo shipper:

- LIBERO PDF Logger Te1-P or Te1-PY
- Temperature Probe
- Stainless steel Bracket
- Mounting fee

Send us an e-mail to swiss@elpro.com indicating your container type and request a quote.



Optional:

elproVIEWER Software
for additional data evaluation of
embedded data in PDF File (free)
or data overlay/join (liable to cost).





Basic Edition:

- View embedded data of all LIBERO PDF Loggers
- Detailed analysis & powerful reports

Professional Edition:

- Assessment functionality incl. specific report
- Overlay of up to 100 curves
- Join up of up to 16 curves



Our Contact Details

ELPRO-BUCHS AG | Langäulistrasse 45
9470 Buchs SG | Switzerland
T +41 81 552 08 08 | www.elpro.com

Additional onsite ELPRO support offices in Benelux and Denmark
Worldwide distributor network